



Product Catalog

2011-10-12 / BG0143-Z



Trademarks

Adobe®

CIDEX®

CIV-Flex™

DICOM®

IQPAC™

Magnum®

Microsoft®

Micro-Touch®

NeoGuard®

Storz®

Welch Allyn®



Table of Contents

Systems	1	Transducer Holder.....	11
Pro Focus 2202 System	1		
Available versions for Pro Focus	1		
Flex Focus 1202 System.....	1		
Available versions for Flex Focus.....	1		
2202 Upgrade Options	3		
Accessories for Pro Focus 2202.....	7	3D Imaging Accessories for	
PCU.....	7	Pro Focus 2202 and Flex Focus 1202	12
Trackball Adjustment Tool.....	7	Magnetic Wheel Mover.....	12
Foot Switch	7	Adapter Ring for 8848.....	12
Hand Rest for 2202	7	Adapter Ring for 8808.....	12
Foot Switch	7	Adapter Ring for 8818/8808e.....	12
Cord Retaining Kit.....	7	Pullback Adapter for 8848	12
Network Isolator.....	7	Pullback Adapter for 8818/8808e.....	12
Transducer Holder Inlay.....	7	Adapter Cable for 1202.....	12
Extension Cable for Front Panel.....	7	Variable Friction Support Arm.....	12
Endo Transducer Holder for Pro Focus.....	8	Foot Switch.....	12
Transducer Holder for Transducer 8838 for use on		Foot Switch.....	12
Pro Focus 2202 UltraView.....	8	3D Connector Box.....	12
Shortcut Stickers for Pro Focus.....	8		
System Cover.....	8	Documentation Accessories	13
Flight Case for 2202.....	8	Printing Paper SONY UPP-110HA	13
Digital B/W Printer Start-Up Kit.....	8	Printing Paper SONY UPP-110HG	13
Printer Shelf.....	8	Color Printing Pack SONY UPC-21S	13
SONY Monitor	8	Digital B/W Printer SONY	
		UP-D897	13
		Digital B/W Printer SONY	
		UP-D897 for Flex Focus 1202	13
		Digital Color Printer SONY	
		UP-D25MD for Pro Focus 2202 UltraView.....	13
		Printer Start-Up Kit for Flex Focus 1202.....	13
Accessories for Flex Focus 1202	9	Power Cables.....	14
Power Pack.....	9		
Battery Charger Station.....	9		
Mobile Keyboard Dock.....	9	Miscellaneous	18
Keyboard Wall Dock.....	9	Scanning Gel	18
Mobile Keyboard Dock.....	9	OB Stickers for keyboard shortcuts.....	18
Keyboard Table Dock.....	9	Breast Ultrasound Biopsy	
Mobile Touch Dock.....	9	Phantom. Model 052	18
Touch Wall Dock	9	Tissue Equivalent Ultrasound Prostate Phantom	
Mobile Touch Dock.....	10	Model 053.....	18
Backpack for Flex Focus.....	10	Biplane Prostate Phantom.....	18
DVD Drive (upgrade).....	10	Endfire Prostate Phantom	18
Remote Control	10		
Wheeled Transport Case	10	Transducers and Transducer	
Endo Transducer Holder	10	Accessories	19
Endo Transducer Holder Mount.....	10	Endosonic Transducer 1850-M.....	19
Wall-Mounted Transducer Holder	10	Endosonic Transducer 1850-M (continued).....	20
Transducer Connector to mount on the wall (dummy)	10	Anorectal Transducers 2050 and 2052	21
Foot Switch	11		
Video Cable (Flex Focus to phono and BNC)	11		



Table of Contents

Mechanical Sector Transducer 8551-R.....	22
Laparoscopic Transducer 8666-RF	22
Convex Array Transducer 8667	23
Linear Array Transducer 8670.....	23
Convex Array Transducer 8802.....	24
Biplane Transducer 8808	25
Endorectal Biplane Transducer 8808e	26
Linear Array Transducer 8809.....	27
Linear Array Transducer 8811.....	27
Linear Array Transducer 8812.....	28
Intraoperative Transducer 8814.....	28
Intraoperative Transducer 8815.....	29
T-Shaped Intraoperative Transducer 8816	29
Biplane and Endfire Transducer 8818	30
Endovaginal Endfire Transducer 8819	31
Convex Array Transducer 8820e.....	31
Vascular Transducer 8822	32
Convex Array Transducer 8823.....	32
I and T Shaped Intraoperative Transducer 8824.....	33
ProART Robotic Transducer 8826.....	33
Convex Array Transducer 8830.....	34
3DART Transducer 8838.....	35
Transducer 8848.....	36
Craniotomy Transducer 8862.....	36
Burr Hole Transducer 8863.....	37
High Frequency Linear Array Transducer 8870.....	37
CIVCO AccuCare Stepper Units (Barzell-Whitmore)	38

Systems

2202



Pro Focus 2202 System

Printer must be ordered separately.

2202-7-UltraView

Available versions for Pro Focus

High-End Color System

2202-2-Green (C)

Midrange Color System

2202-9-OR

Pro Focus UltraView OR

1202



Flex Focus 1202 System

Keyboard Dock, printer etc. must be ordered separately.

1202-3

Available versions for Flex Focus

Flex Focus 700

Dedicated to surgery.

1202-1

Flex Focus 400

Color system with IQPAC™.

1202-4

Flex Focus 400 Anesthesia

Dedicated to regional anesthesia.

1202-2

Flex Focus 200

Flex Focus with basic functionality. B&W system. Upgradable to color.

1202-7

Flex Focus 200 Surgery

Basic surgical system with limited transducer range.

Systems

Available versions of Pro Focus 2202		
Part no. Pre-order	Feature	UltraView 2202-7
	IQPAC™	Included
	Color Flow Mapping	Included
	Power Doppler	Included
	Directional Power Doppler	Included
	Pulsed Wave Doppler	Included
	Continuous Wave Doppler	Included
UA2214	Tissue Harmonic Imaging	Included
UA2215	Contrast Imaging (License) ¹	Optional
UA2206	DICOM Option (License)	Optional
UA1205	Picture-in-Picture	Optional
UA2238	HistoScanning™ Interface Option for 2202 UltraView ²	Optional
UA2238-K	HistoScanning™ Interface Upgrade Kit for 2202 UltraView ²	Optional
UA2235	HistoScanning™ Interface Option for 2202 Yellow (CX) ²	NA
UA2235-K	HistoScanning™ Interface Upgrade Kit for 2202 Yellow (CX) ²	NA
UA2236	HistoScanning™ Interface Option for 2202 Green (C) ²	NA
UA2236-K	HistoScanning™ Interface Upgrade Kit for 2202 Green (C) ²	NA
UA2237-K	HistoScanning™ Interface Upgrade Kit for 2202 Yellow (CX) and Green (C) ²	NA
UA2226	Varian Option License	Optional
UA2226-K	Varian Option Upgrade Kit	Optional
UA2207	3D Freehand (License)	Optional
UA2222	3D Professional (includes 3D Mover Module and Freehand)	Optional
UA2222-K	3D Mover Module Kit (from Freehand to Professional. Must be ordered with UA2207)	Optional
UA2223	360° Mechanical Input Module	Optional
UA2223-K	360° Mechanical Input Module Kit (Upgrade)	Optional
UA2251	Flat Screen for 2202	NA
UA2251-K	Flat Screen Kit for 2202	NA
UA2252	Flat Screen on Arm for 2202	NA
UA2252-K	Flat Screen on Arm Kit for 2202	NA
	Height-Adjustable Keyboard Panel and Monitor (Simultaneous Movement)	Included
	Height-Adjustable Keyboard Panel and Monitor (Independent Movement)	NA
UA2220	3rd Linear and Convex Array (Electronic) Transducer Connector Active	Included

Systems

Available versions of Pro Focus 2202		
Part no. Pre-order	Feature	UltraView 2202-7
	Support for transducers: 8667, 8670, 8802, 8808, 8809, 8811, 8818, 8819, 8820e, 8823, 8827, 8848	Included
	Support for transducers: 8666-RF, 8808e, 8812, 8814, 8815, 8816, 8822, 8824, 8826, 8838, 8862, 8863, 8870	Included
	Support for Single Element (Mechanical) transducers: 1850, 2050, 2052, 8551 ³	Optional
UA2410	Palm Control Unit (PCU). See page 7 for more details.	Optional
UA4100	Digital B/W Printer Start-Up Kit. See page 8 for more details.	Optional

1 In the USA, contrast-enhanced ultrasound has not been market cleared by the FDA, with the exception of only select cardiac applications.

2 HistoScanning has not been market cleared by the FDA for sale in the USA.

3 Transducer Type 8551 is not supported on UltraView.

2202 Upgrade Options

UA2233-K

B/W to Color Upgrade Kit[†]

UA2229-K

Transducer Support Upgrade Kit
Required to support 8818 on systems from before March 2006 and to support 8848 on systems from before January 2007. Contact Service for further details.

[†] For systems purchased after February 2007.

Systems

Options for Flex Focus 1202						
Part no. Pre-order	• Included ° Optional Feature	Flex Focus 200	Flex Focus 200 Surgery	Flex Focus 400	Flex Focus 400 Anesthe- sia	Flex Focus 700
UA1210	Mobile Keyboard Dock	°	°	°	°	°
UA1214	Mobile Keyboard Dock with battery (includes 2 rechargeable batteries)	°	°	°	°	°
UA1211	Keyboard Wall Dock	°	°	°	°	°
UA1283	Keyboard Table Dock	°	°	°	°	°
UA1810	Mobile Touch Dock	°	°	°	°	°
UA1811	Touch Wall Dock (no DVD [RW])	°	°	°	°	°
UA1814	Mobile Touch Dock with battery (includes 2 rechargeable batteries)	°	°	°	°	°
UA1222	3D Professional (License) - includes 3D Freehand	NA	NA	°	NA	°
UA1207	3D Freehand (License)	NA	NA	°	°	°
UA1206	DICOM (License)	°	°	°	°	°
UA1226	Varian Option (License)	NA	NA	°	NA	°
UA1203	Color Doppler (Color Flow Mapping)	1	°	•	•	•
UA1298	Pulsed Wave Dopple (PW)	1	°	•	•	•
UA1205	Picture-in-Picture	NA	NA	NA	NA	°
	IQPAC™	NA	NA	•	•	•
	Clip Store	NA	NA	•	•	•
	Video Out	NA	•	•	•	•
	Basic Transducer Support: 8667, 8670, 8808, 8808e, 8811, 8819, 8820e, 8823, 8827, 8830	•	•	•	NA	•
	Enhanced Transducer Support: 1850, 2050, 2052, 8802, 8809, 8818, 8822, 8838, 8848, 8862, 8863, 8870	NA	NA	•	NA	•
	Anesthesia Transducer Support: 8670, 8809, 8811, 8816, 8823, 8827, 8830, 8870	NA	NA	NA	•	NA
UA1258	Basic Surgical Transducer Support: 8666-RF, 8815, 8816	NA	•	NA	NA	NA
	Complete Surgical Transducer Support (all transducers)	NA	NA	NA	NA	•
UA1237	Remote Control for Flex Focus 700 (includes Mini Bluetooth Adapter)	NA	NA	NA	NA	°
UA1294	Endo Transducer Holder	NA	°	°	NA	°
UA1295	Endo Transducer Holder Mount (for UA1294 with 1850 or 2052)	NA	°	°	NA	°
UA1204	Wall-Mounted Transducer Holder	°	°	°	°	°

Systems


Options for Flex Focus 1202

Part no. Pre-order	• Included ° Optional Feature	Flex Focus 200	Flex Focus 200 Surgery	Flex Focus 400	Flex Focus 400 Anesthe- sia	Flex Focus 700
UA1296	General Transducer Holder (requires UA1810 or UA1814)	°	°	°	°	°
UA1223	Power Cable Holder. Included with UA1214 and UA1814	°	°	°	°	°
UA1220	Backpack to hold printer. Requires Mobile Dock UA1210	°	°	°	°	°
UA1202	Foot Switch (USB)	°	°	°	°	°
UA1284-K	DVD RW Drive Kit	°	°	°	°	NA
UA1247	Battery Charger Station ¹	°	°	°	°	°
UA1225	Extra Battery for UA 1214 and UA1814 (rechargable 14.4 V lithium ion battery). Must be ordered in sets of 2 or 4 batteries.	°	°	°	°	°
QA0228	Trackball Adjustment Tool	•	•	•	•	NA
UA4102	Printer Start-Up Kit. Requires UA1220, UA1214 or UA1814	°	°	°	°	°
AO0491	Video Cable (Flex Focus to phono and BNC)	°	°	°	°	NA
UA0964	Flight Case - for scanner unit with Mobile Keyboard Dock	°	°	°	°	°
UA1285	Wheeled Transport Case - suitcase on wheels for scanner unit without dock, plus two transducers.	°	°	°	°	°

¹ Only available in some countries.

flex Focus

Flexibility concept



Accessories for Pro Focus 2202

UA2410



PCU

Palm Control Unit for 2202. Allows the user to operate all major functions remotely. The Palm Control Unit may be disinfected and sterilized. Ideal for the surgical environment.

UA1271

Adapter for PCU. Needed to connect UA2410 to 2202.

UA0076

Sterile cover, 14 x 147 cm, latex-free

UA1202



Foot Switch

For Pro Focus 2202 Ultraview and Flex Focus 1202.

QA0228



Trackball Adjustment Tool

Tool to adjust the trackball tension on 2202 and 1202.

Included with the system.

UA1215



Foot Switch

2-pedal high-quality foot switch. Both pedals are user-configurable when connected to 2202 Yellow, Green or Blue.

DS0967



Hand Rest for 2202

The removable hand rest can be used to make working with the keyboard panel more comfortable.

Cables and Connectors

UA2204

Cord Retaining Kit

DP0925

Network Isolator

DZ9756



Transducer Holder Inlay

Inlay for the transducer holders on the system. This allows very small transducers, such as the 8809, to be securely stored on the system. Only if extras are needed. 1 is included as standard with the transducer.

AO0421



Extension Cable for Front Panel

Unique BK Medical feature. Using the extension cable, the system can be operated from up to 3 meters away.

Accessories for Pro Focus 2202

UA1302



Endo Transducer Holder for Pro Focus


Holds all large endo transducers and the ECRM-Mover in a stable position when not in use.

SONY Monitor

SONY LMD2140MD medical grade LCD monitor.



UA1347



Transducer Holder for Transducer 8838 for use on Pro Focus 2202 UltraView.

SC1624

Shortcut Stickers for Pro Focus

UA0955

System Cover

UA0959

Flight Case for 2202

UA4100

Digital B/W Printer Start-Up Kit

Must be ordered separately. Consists of a SONY UP-D897 printer, a USB cable, a power cable and 10 rolls of printer paper.

UA2240

Printer Shelf

Accessories for Flex Focus 1202

Flex Focus Scanner Units

1201-1, 2, 4 or 7



1202-1 (Flex Focus 400) is a color scanner unit. The 1202-2 (Flex Focus 200) is a Flex Focus scanner unit, with basic functionality and 1202-4 (Flex Focus 400 Anesthesia) is dedicated to regional anesthesia. A Flex Focus 200 system or 400 Anesthesia can only be upgraded by replacing the scanner unit with a Flex Focus 400 scanner unit. To get a complete Flex Focus system, you need to order both a scanner unit and a dock. 1202-7 (Flex Focus 200 Surgery) is an affordable surgical system with basic capabilities and limited surgical transducer support.

1201-3



1202-3 (Flex Focus 700) is a dedicated surgical system with touch screen. You need to order both a scanner unit and a dock.

Power Pack



The Power Pack enables you to operate the system while it is not connected to an external electrical supply. It is located in the bottom drawer of the backpack. Only with Mobile Battery Keyboard Dock with battery UA1214 or Mobile Touch Dock with Battery UA1814.

UA1247




Battery Charger Station

Like the Power Pack, but designed to stand on a table or hang on a wall.


Docks

UA1210



Mobile Keyboard Dock

UA1211




Keyboard Wall Dock

UA1214



Mobile Keyboard Dock
With battery

UA1283



Keyboard Table Dock

UA1810



Mobile Touch Dock
No DVD (RW) option

UA1811



Touch Wall Dock
No DVD (RW) option

Accessories for Flex Focus 1202

Docks

UA1814



Mobile Touch Dock
With battery

UA1285



Wheeled Transport Case

Suitcase on wheels for 1202. Holds:

- Scanner unit without dock
- Two transducers

UA1220



Backpack for Flex Focus

Transducer connector parking on top. Holds either two drawers, one drawer and a printer or the Battery Pack.

UA1294



Endo Transducer Holder

For mount on the array connector.

UA1284-K



DVD Drive (upgrade)

Requires Keyboard Dock

UA1295



Endo Transducer Holder Mount

Dummy array connector for mechanical transducers without array connectors.

If you want to use Endo Transducer Holder UA1294 for 1850, 2050 or 2052, UA1295 is needed.

UA1237



Remote Control

For Flex Focus 700. Keys and mouse function to operate the system. Runs on two AA batteries. Includes Mini Bluetooth Adapter.

UA1204



Wall-Mounted Transducer Holder Transducer Connector to mount on the wall (dummy)

Together with UA1204 this allows transducer storage on the wall, for instance near the sink.

Accessories for Flex Focus 1202

UA1308

Transducer Holder

For Flex Focus 1202



UA1202

Foot Switch

For Pro Focus 2202 Ultraview and Flex Focus 1202



AO0491

Video Cable (Flex Focus to phono and BNC)

Allows you to connect Flex Focus to video devices with phono and BNC sockets. 2.5 m long.



3D Imaging Accessories for Pro Focus 2202 and Flex Focus 1202

UA0513



Magnetic Wheel Mover

Available for use with transducer types 8658, 8808, 8808e, 8818 and 8848

The mover is designed to rotate rectal ultrasound transducers in order to acquire a dataset that can be used to construct a 3D volume. It can be used as both a rotational mover and a pull-back mover.

UA0510

Adapter Ring for 8848

For use with UA0513 as a rotational mover.

UA0511

Adapter Ring for 8808

For use with UA0513 as a rotational mover.

UA0512

Adapter Ring for 8818/8808e

For use with UA0513 as a rotational mover.

UA0515

Pullback Adapter for 8848

For use with UA0513 as a pullback mover.

UA0516

Pullback Adapter for 8818/8808e

For use with UA0513 as a pullback mover.

EL4036



Adapter Cable for 1202

For use with UA0513 as a pullback mover.

UA1215



Foot Switch

2-Pedal High-Quality Foot Switch. Both pedals are user-configurable when connected to 2202.

UA1202



Foot Switch

For Pro Focus 2202 Ultraview and Flex Focus 1202

UA0670

3D Connector Box

Required for 3D connection to the following systems: 6503, 6503-2, 6503-4, 2202 with UA2202*.

*2202 is currently delivered with UA2223 which does not require a 3D Connector Box.

UA0553



Variable Friction Support Arm

Variable Friction Arm for UA0513 (Magnetic Wheel Mover)

Documentation Accessories

EQ4042



Printing Paper SONY UPP-110HA

High density printing paper.
Roll size: 110 mm x 18 m.

For older printers.
Printer is no longer available.

EQ4104



Digital B/W Printer SONY UP-D897 for Flex Focus 1202

Digital Monochrome printer with a USB interface. The UP-D897 is a feature packed printer to allow for the greatest flexibility in black and white prints.

Part of Printer Start-Up Kit UA4102
Paper: EQ4072

EQ4072



Printing Paper SONY UPP-110HG

High glossy printing paper.
Roll size: 110 mm x 20 m.
Print quantity: Approx. 217 prints per roll. 10 rolls of paper included with the printer.

EQ4105



Digital Color Printer SONY UP-D25MD for Pro Focus 2202 UltraView

Digital color dye sublimation printer with a USB interface. The UP-D25MD features a slim, compact design that is ideal for medical carts.

Paper: EQ4093 and EQ4094

EQ4093



Color Printing Pack SONY UPC-21S

240 sheets of print paper.
3 rolls of ink ribbon.
Paper size: 100 x 90 mm.

UA4102

Printer Start-Up Kit for Flex Focus 1202

Must be orderede separately. Consists of a SONY UP-D897 printer (EQ4104), a USB cable, a power cable and 10 rolls of printer paper.

EQ4100



Digital B/W Printer SONY UP-D897


Digital Monochrome printer with a USB interface. The UP-D897 is a feature packed printer to allow for the greatest flexibility in black and white prints.

Part of Printer Start-Up Kit UA4100
Paper: EQ4072

Power Cables

AN0084

2202 and 1202




AN0089

2202 and 1202




AN0082

2202 and 1202




AN0081

2202 and 1202




AN0083

2202 and 1202




AN0080

2202 and 1202



AN0090


2202 and 1202



(Reversed phase and null)

AN0086


2202 and 1202



Round pin


AN0085

2202 and 1202



AN0087


2202 and 1202



Flat pin

AN0088

2202 and 1202



Power Cables

COUNTRY	SIGNAL TYPE	COUNTRY CODE	2202/1202	VOLTAGE
ALBANIA	PAL B/G	915	AN0080	230 V
ALGERIA	PAL B/G	903	AN0080	230 V
ARGENTINA	PAL N	921	AN0090	230 V
AUSTRALIA	PAL B/G	918	AN0083	230 V
AUSTRIA	PAL B/G	904	AN0080	230 V
AZERBAIJAN	PAL B/G	915	AN0080	230 V
BAHRAIN	PAL B/G	900	AN0082	230 V
BANGLADESH	PAL B	900	AN0082	230 V
BELARUS	SECAM D/K	915	AN0080	230 V
BELGIUM (LUXEMBOURG)	PAL B/H	901	AN0080	230 V
BOLIVIA	NTSC M	913	AN0081	230 V
BOSNIA/HERZEGOVINA/ SERBIA	PAL B/H	915	AN0080	230 V
BRAZIL	PAL M	914	AN0081	115 V
BULGARIA	SECAM D/K	992	AN0080	230 V
CAMBODIA	"PAL B/G, NTSC M"	900	AN0082	230 V
CHILE	NTSC M	923	AN0085	230 V
CHINA (PEOPLE'S REPUB- LIC)	PAL D	906	AN0083	230 V
COLOMBIA	NTSC M	913	AN0081	115 V
COSTA RICA	NTSC M	913	AN0081	115 V
CROATIA	PAL B/H	915	AN0080	230 V
CUBA	NTSC M	913	AN0081	115 V
CYPRUS	PAL B/G	983	AN0082	230 V
CZECH REPUBLIC	"PAL B/G, SECAM D/K"	982	AN0080	230 V
DENMARK	PAL B/G	974	AN0087/AN0086	230 v
ECUADOR	NTSC M	913	AN0081	115 V
EGYPT	"PAL B/G, SECAM B/G"	900	AN0082	230 V
EL SALVADOR	NTSC M	913	AN0081	115 V
ESTONIA	"PAL B, SECAM D/K"	984	AN0080	230 V
FINLAND	PAL B/G	908	AN0080	230 V
FRANCE	SECAM L	903	AN0080	230 V
GEORGIA	SECAM D/K	915	AN0080	230 V
GERMANY	PAL B/G	904	AN0080	230 V
GHANA	PAL B/G	900	AN0082	230 V
GREECE	SECAM B/G	910	AN0080	230 V
GREENLAND	PAL B	974	AN0087/AN0086	230 V
GUATEMALA	NTSC M	913	AN0081	115 V
HONDURAS	NTSC M	913	AN0081	115 V
HONG KONG	PAL I	924	AN0082	230 V
HUNGARY	PAL K/K	987	AN0080	230 V
ICELAND	PAL B/G	915	AN0080	230 V
INDIA	PAL B	917	AN0084	230 V
INDONESIA	PAL B	915	AN0080	230 V
IRAN	SECAM B/G	915	AN0080	230 V

Power Cables

COUNTRY	SIGNAL TYPE	COUNTRY CODE	2202/1202	VOLTAGE
IRAQ	SECAM B/G	900	AN0082	230 V
ISRAEL	PAL B/G	919	AN0089	230 V
ITALY	PAL B/G	902	AN0085	230 V
JAPAN	NTSC M	905	AN0081	115 V
JORDAN	PAL B/G	900	AN0082	230 V
KAZAKHSTAN	SECAM D/K	915	AN0080	230 V
KENYA	PAL B/G	900	AN0082	230 V
KUWAIT	PAL B/G	900	AN0082	230 V
LATVIA	"PAL B/G, SECAM D/K"	985	AN0080	230 V
LEBANON	PAL B/G	900	AN0082	230 V
LITHUANIA	"PAL B/G, SECAM D/K"	986	AN0080	230 V
MALAYSIA	PAL B	900	AN0082	230 V
MALTA	PAL B	988	AN0082	230 V
MEXICO	NTSC M	913	AN0081	115 V
MOROCCO	SECAM B	903	AN0080	230 V
NETHERLANDS	PAL B/G	901	AN0080	230 V
NEW ZEALAND	PAL B/G	918	AN0083	230 V
NIGERIA	PAL B/G	900	AN0082	230 V
NORWAY	PAL B/G	976	AN0080	230 V
OMAN	PAL B/G	900	AN0082	230 V
PAKISTAN	PAL B	900	AN0082	230 V
PERU	NTSC M	907	AN0080	230 V
PHILIPPINES	NTSC M	916	AN0081	230 V
POLAND	PAL B/G	981	AN0080	230 V
PORTUGAL	PAL B/G	911	AN0080	230 V
QATAR	PAL B	900	AN0082	230 V
ROMANIA	PAL D/G	991	AN0080	230 V
RUSSIA	SECAM D/K	977	AN0080	230 V
SAUDI ARABIA	"SECAM B/G, PAL B"	900	AN0082	230 V
SINGAPORE	PAL B/G	900	AN0082	230 V
SLOVAKIA	"PAL B, SECAM D/K"	990	AN0080	230 V
SLOVENIA	PAL B/H	989	AN0080	230 V
SOUTH AFRICA	PAL I	917	AN0084	230 V
SOUTH KOREA	NTSC M	922	AN0080	230 V
SPAIN	PAL B/G	907	AN0080	230 V
SWEDEN	PAL B/G	975	AN0080	230 V
SWITZERLAND	PAL B/G	920	AN0088	230 V
SYRIA	"SECAM B, PAL G"	915	AN0080	230 V
TAIWAN	NTSC	905	AN0081	115 V
THAILAND	PAL B/M	916	AN0081	230 V
TUNISIA	SECAM B/G	903	AN0080	230 V
TURKEY	PAL B	915	AN0080	230 V
UKRAINE	SECAM D/K	915	AN0080	230 V
URUGUAY	PAL N	915	AN0080	230 V
UNITED ARAB EMIRATES	PAL B/G	900	AN0082	230 V

Power Cables

COUNTRY	SIGNAL TYPE	COUNTRY CODE	2202/1202	VOLTAGE
UNITED STATES (CANADA)	NTSC M	909	AN0081	115 V
UNITED KINGDOM	PAL I	900	AN0082	230 V
VENEZUELA	NTSC M	913	AN0081	115 V
VIETNAM	"NTSC M, SECAM D"	915	AN0080	230 V
YEMEN	PAL I	915	AN0082	230 V

Miscellaneous

WQ0739



Scanning Gel

Electronic coupling medium for ultrasound imaging.
Refillable dispenser 290 ml.


Refill container 5 liters:
UA0610

WQ0970

Biplane Prostate Phantom

Training and simulation device.
Organs and tissues are correctly adjusted to mimic exact Ultrasound echogenic seen by surgeon during OR procedures.

SC1474



OB Stickers for keyboard shortcuts

Extremely useful for users who use the programmable shortcuts in the obstetrics package to program shortcuts for advanced measurements. The stickers are placed on the shortcut keys you define.

WQ0971

Endfire Prostate Phantom

Training and simulation device.
Organs and tissues are correctly adjusted to mimic exact Ultrasound echogenic seen by surgeon during OR procedures.

WQ0975



Breast Ultrasound Biopsy Phantom. Model 052

Training device for ultrasound guided biopsy procedures.
The phantom mimics the ultrasonic characteristics of tissues found in an average human breast.


WQ0974



Tissue Equivalent Ultrasound Prostate Phantom. Model 053

Training device for imaging the prostate with a rectal transducer.


Transducers and Transducer Accessories



Endosonic Transducer 1850-M

Systems: 2202 with Single Element Module UA2223/UA2223-K, 1202-1 and 1202-3

Frequency Range: 10-5 MHz Includes as standard: 1 transducer head, 1 tube and 1 hood. Must be specified when ordering. Type 1850-M without accessories: Order 1850-X or 1850-R (replacement)



Rectal Transducer

Type 6004

For transrectal examination.

10 MHz, 90°



UA0543

17 mm Rectal Tube
Perineal Puncture Attachment
Hood for UA0878
Hood for UA0878, plastic

UA0878
UA0650
WA0543
EA4020

Incl. lubricant and Rubber Sheaths UA0799.
Bore diameter 2.1 mm.
Autoclavable.



EA4020

Standoff Sheaths (box of 100)

UA0799

Non-sterile for use with Rectal Tube UA0878.
1 box included with UA0878.



UA0681

Inner Cylindrical Ring for UA0878
Outer Cylindrical Ring for UA0878
O-Rings for WA0543/EA4020

YA0514
YM1615
YJ0621

Sterile, Single-Use Rectosigmoidoscope

UA0681

With obturator and grip.




UA0682

Light Source Cable Adapter for UA0681

UA0682

Reusable. For Storz and Wolf Light Sources.



Urethral and Rectal Transducer Type 6005

For transurethral and pediatric endoanal examination.

10, 7.0 and 5.0 MHz, 90°

Tube for Type 6005, 10 mm
Extra hoods for WA0500, convex
Extra hoods for WA0500, slim

WA0500
WA0499
WA0597

Incl. 2 hoods and Lubricant QS0011.

Extra hoods for WA0637

WS3416


Storz Sheath
Adapter for WQ0410

WQ0410
UA0652

Transducers and Transducer Accessories



Endosonic Transducer 1850-M (continued)



Type 8539

Rectal Transducer

10, 7.0 and 5.0 MHz, 90°

Rectal Transducer
5.0, 7.0 and 10 MHz,

Type 8539

For transrectal examination.
Same accessories as for Type 6004.



Perineal Puncture Attachment
Standoff Sheaths (box of 100)
Lubricant

UA0650
UA0799
QS0011

Bore diameter 2.1 mm.
Non-sterile for use with Rectal Tube UA0878.

O-Ring for sealing between 1850-M and
Rectal Tubes
O-Rings for crystals

YJ0632
YJ0619

(proximal end of 1850-M)



Rasha Applicator

EA4006

Latex Covers

UA0069

24 pcs., Sterile, 3.2 cm tapered to 1.9 x 32 cm.

Transducers and Transducer Accessories

Anorectal Transducers 2050 and 2052

Systems: 2202 with Single Element Module UA2223/UA2223-K, 1202-1 and 1202-3
Frequency Range: 16-6 MHz



UA0673



UA0681




UA0679

Carrying Case	KE4322	Included.
Water Standoff System:		
Water Standoff Collar	UA0671	
Rectosigmoidoscope with obturator	UA0672	For Welch Allyn® Light Source.
Rectosigmoidoscope with obturator	UA0673	For Storz® Light Source.
Rectosigmoidoscope. Sterile, single-use*	UA0681	With obturator and grip.
Non-Sterile Water Standoff Condoms, latex	UA0037	24 covers 2.6 x 30 cm.
Sterile Latex-Free Water Standoff Condoms	UA0068	24 pcs.
Light Source Cable Adapter for UA0681*	UA0682	Reusable. For Storz and Wolf Light Sources.
Spare O-Rings for use with UA0671	UA0674	Pack of 10 Silicone O-Rings.
Glycerin-Based Lubricant	UA0675	250 ml for UA0037.
Plastic Extension Tube, 30 cm	UA0676	(Included with UA0671 - can also be purchased separately).
2-Way Stopcock	UA0677	(Included with UA0671 - can also be purchased separately).
Plastic Syringe	UA0678	(Included with UA0671 - can also be purchased separately).
CIDEX® Instrument Tray System	UA0061	
3D Imaging:		
3D Connector Box	UA0670	Required for 3D connection to the following systems: 6503, 6503-2, 6503-4, 2202 with UA2202**.
Variable Friction Support Arm	UA0553	
Transducer Support Collar	UA0679	For use with Variable Friction Support Arm UA0553 or any other clamp.
Sterile Latex Transducer Covers	UA0005	24 pcs. sterile transducer cover 2.6 x 20 cm.

*Not released for sale in Canada

**2202 is delivered with UA2223 which does not require a 3D Connector Box.

Transducers and Transducer Accessories



Mechanical Sector Transducer 8551-R

Systems: 2101 Premium with Single Element Module UA1255 or UA1255-K, 1101 with an integrated Single Element Module and 2202 Yellow, Green and Blue with Single Element Module, Frequency Range: 10 - 6 MHz Only sold as a replacement transducer

Carrying Case	KE4300	Included.
Insert for Carrying Case	DS0524	Included.
Extended Needle Guide for transrectal puncture	UC5299-E	114.5 mm long.
Single-Use Cleaning Brush	QZ0040	For needle guide. Included.
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
Latex Covers	UA0006	24 pcs., sterile covers 3.5 x 20 cm.
Latex Covers	UA0069	24 pcs., sterile, 3.2 cm tapered to 1.9 x 32 cm.
Water Stand-Off Rubber Sheath	UA0799	1 Pack of 100 included. Non-sterile.
Water Stand-Off O-Ring	QA0185	1 Pack of 100 included. Non-sterile.



Laparoscopic Transducer 8666-RF

Systems: 2202, 1202-3 and 1202-7
Frequency Range: 10-5 MHz



UA1166



UA1249

Carrying Case	KE4320	Included.
Biopsy Attachment	UA1249	Needed for UA0039.
Histological Biopsy Needle Kit	UA0039	19 gauge needles needed. Box of 5 pcs.
No-Flex™ Attachment	UA1166	For customers who prefer to image with a rigid laparoscopic transducer.
Single-Use Cleaning Brush, thin	QZ0040	For needle guide.
NeoGuard® Sterile Covers	UA0067	12 pcs, telescopically folded, latex-free.
Leakage Testing Kit	UA1404	

Transducers and Transducer Accessories



Convex Array Transducer 8667

Systems: 2202, 1202-1, 1202-2, 1202-3 and 1202-7
Frequency Range: 10-5 MHz

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1256	Bore diameter of 1.5 mm and angle of insertion 0° to the long axis of the transducer.
Single-Use, Sterile Needle Guide	UA1256-U*	Bore diameter of 1.3 mm and angle of insertion 0° to the long axis of the transducer. 8 pcs.
Leakage Testing Kit	UA1404	
Latex Covers	UA0006	24 pcs. sterile covers 3.5 x 20 cm.
CIV-Flex™ Latex-Free, Sterile Covers	UA0008	24 pcs 11.97 cm narrowing down to 4.57 cm.

* UA1256-U not released for sale in the USA and Canada



Linear Array Transducer 8670

Systems: 2202 and 1202
Frequency Range: 12-4 MHz

Carrying Case	KE4300	Must be ordered separately.	
UA1246	Needle Guide (for 8570 with 1101 and for 8670 with 1402 only)	UA1246	Bore diameter 0.7-2.4 mm and angle of insertion of 45° to the image axis. Requires UA0028.
	Needle Guide (For 8670 with 2202 and 1202 only)	UA1269	Bore diameter 0.6-2.4 mm and angle of insertion of 30°, 45° and 60° to the image axis.
	Needle Guide Insert for UA1246	UA0028	Box with 20 single-use inserts.
UA1269	Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
	CIV-Flex™ Sterile, Latex-Free Covers	UA0009	24 pcs. 8.89 x 60 cm.
	CIV-Flex™ Sterile Covers, latex-free NeoGuard® Sterile Covers	UA0071 UA0074	24 pcs., 8.9 x 91.5 cm, telescopically folded. 12 pcs., latex-free, 20.3 x 244 cm with SurgiTip 6 cm, telescopically folded.

Transducers and Transducer Accessories



Convex Array Transducer 8802

Systems: 2202, 1202-1 and 1202-3
Frequency Range: 6-3 MHz



UA1250

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.
Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009	24 pcs. 8.89 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0071	24 pcs., 8.9 x 91.5 cm, telescopically folded.
CIV-Flex™ Cover, latex-free	UA0076	24 pcs., 14 x 147 cm, telescopically folded.

Transducers and Transducer Accessories



Biplane Transducer 8808

Systems: 2202, 1202-1, 1202-2, 1202-3 and 1202-7
Frequency Range: 10-5 MHz

Carrying Case	KE4300	Must be ordered separately.
Dummy Channel Bracket	UA1272	Included.

Transrectal Biopsy Guides

Reusable Biopsy Guide	UA1257-R17 UA1257 UC5296	For 17-gauge needles. Includes: Biopsy Channel Bracket 17-gauge needle guide
Extended, Reusable Biopsy Guides	UA1257R17E UA1257 UC5296-E	For 17-gauge needles. Long needle guide. Includes: Biopsy Channel Bracket 17-gauge long needle guide
Reusable Biopsy Guide	UA1257-R14 UA1257 UC5298	For 14-gauge needles. Includes: Biopsy Channel Bracket 14-gauge needle guide



UA1257-S14



UA1257S17E

Single-Use Biopsy Guides	UA1257-S14	For 14-gauge needles. Sterile packed. 18 pcs.
Extended, Single-Use Biopsy Guides	UA1257S17E	For 17-gauge needles. Sterile-packed. 18 pcs. Includes: Biopsy Channel Bracket 17-gauge needle guide



UA1232

Transperineal Puncture Attachment	UA1232	Bore diameter of 2.2 mm and angle of insertion 0° to the long axis of the transducer.
-----------------------------------	--------	---

Single-Use Cleaning Brush, thin	QZ0040	For needle guide.
Magnetic Wheel Mover for 3D Imaging Mover Adapter Ring for UA0513 Variable Friction Support Arm	UA0513 UA0511 UA0553	For Magnetic Wheel Mover.
Rubber Water Standoff Sheaths Water Standoff O-Rings	UA0799 UA0210	Non-sterile. Pack of 100. Non-sterile. Pack of 100.
Leakage Testing Kit CIV-Flex™ Latex-Free, Sterile Covers	UA1404 UA0008	24 pcs. 11.97 cm narrowing down to 4.57 cm.
Latex Covers Latex Covers	UA0032 UA0069	24 pcs., non-sterile. 24 pcs., sterile, 3.2 cm tapered to 1.9 x 32 cm.

Transducers and Transducer Accessories



Endorectal Biplane Transducer 8808e

Systems: 2202 , 1202-1, 1202-2, 1202-3 and 1202-7

Frequency Range: 10-5 MHz



UA 1324

Dummy Channel Bracket	UA1325	Included.
Carrying Case	KE4300	Must be ordered separately.
Simultaneous Biplane Biopsy Guides	UA1326	Reusable. Bore diameter 1.6 mm; angle of insertion 19° to the long axis of the transducer.
Single-Use Biplane Biopsy Guides	UA1322-S	Sterile packed. 18 pcs. Bore diameter of 1.6 mm; angle of insertion 19° to the long axis of the transducer.
Single-Use Biplane Biopsy Guides	UA1322-S14	Sterile packed. 18 pcs. Bore diameter of 2.1 mm; angle of insertion 19° to the long axis of the transducer.
Single-Use Cleaning Brush, thin	QZ0040	For needle guide.
Perineal Puncture Attachment	UA1324	Bore diameter 2.1 mm; angle of insertion 0° to the long axis of the transducer.
Civco Multi-Purpose Workstation	UA2020	For transperineal biopsies.
Magnetic Wheel Mover for 3D Imaging	UA0513	
Mover Adapter Ring for UA0513	UA0512	
Variable Friction Support Arm	UA0553	For Magnetic Wheel Mover.
Pullback Adapter	UA0516	For Magnetic Wheel Mover.
Leakage Testing Kit	UA1404	
CIV-Flex™ Latex-Free, Sterile Covers	UA0008	24 pcs. 11.97 cm narrowing down to 4.57 cm.
Latex Covers	UA0032	24 pcs., non-sterile.
Latex Covers	UA0069	24 pcs., sterile, 3.2 cm tapered to 1.9 x 32 cm.

Transducers and Transducer Accessories



Linear Array Transducer 8809

Systems: 2202 , 1202-1, 1202-3 and 1202-4
Frequency Range: 15-6 MHz



UA1262

Carrying Case	KE4300	Must be ordered separately.
Transducer Holder	DZ9756	Included.
Needle Guide	UA1262	Can be used with needles up to 1.2 mm/18 gauge.
Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004*	24 pcs. 2.5 m.
CIV-Flex™ Latex-Free, Sterile Covers	UA0008*	24 pcs. 11.97 cm narrowing down to 4.57 cm.
NeoGuard® Covers, latex-free	UA0075	12 pcs., 15.2 x 244 cm plus form-fitted NeoGuard® tip, telescopically folded.
* Must be used together		



Linear Array Transducer 8811

Systems: 2202 and 1202
Frequency Range: 12-5 MHz



UA1239

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1239	Bore diameter up to 2.1 mm and angle of insertion 30°, 45° and 60° to the image axis.
Leakage Testing Kit	UA1404	
CIV-Flex™ Covers	UA0007	24 pcs. sterile covers 13.97 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0070	24 pcs., 14 x 91 cm, telescopically folded.
NeoGuard® Sterile Covers	UA0074	12 pcs., latex-free, 20.3 x 244 cm with SurgiTip 6 cm, telescopically folded.

Transducers and Transducer Accessories



Linear Array Transducer 8812

Systems: 2202
Frequency Range: 8-3.75 MHz



UA1246

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1246	Bore diameter of 0.7-2.7 mm and angle of insertion Of 45° to the image axis. Requires UA0028.
Single-Use Needle Cartridges	UA0028	Box of 20. Single-use.
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
CIV-Flex™ sterile, Latex-Free Covers	UA0009	24 pcs. 8.89 x 60 cm.
CIV-Flex™ sterile Covers, latex-free	UA0071	24 pcs., 8.9 x 91.5 cm, telescopically folded.
NeoGuard® Sterile Covers	UA0074	12 pcs., latex-free, 20.3 x 244 cm with SurgiTip 6 cm, telescopically folded.



Intraoperative Transducer 8814

Systems: 2202, 1202-3
Frequency Range: 10-5 MHz



UA1335

Transducer Holder	DZ9756	Included.
Instructions for use	BB1943	Included.
Carrying Case	KE4300	Must be ordered separately.
Single-Use Needle Guide	UA1335	Box of 9.
Leakage Testing Kit	UA1404	
Sterile cable Covers, latex-free	UA0004*	24 pcs. 2.5 m. Used with UA0009.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009*	24 pcs. 8.89 x 60 cm. Used with UA0004.
NeoGuard® Sterile Covers	UA0074	12 pcs., latex-free, 20.3 x 244 cm with SurgiTip. 6 cm, telescopically folded.

*Must be used together

Transducers and Transducer Accessories



Intraoperative Transducer 8815

Systems: 2202, 1202-3 and 1202-7
Frequency Range: 10-5 MHz



UA1336

Transducer Holder	DZ9756	Included.
Carrying Case	KE4300	Must be ordered separately.
Single-Use Needle Guide	UA1336	Includes: 3-angle insert for 10-20 gauge needles (box of 9), Free-angle insert for 10-20 gauge needles (box of 9).
Leakage Testing Kit for 8815	UA1404	
Sterile Cable Covers, latex-free	UA0004*	24 pcs. 2.5 m.
CIV-Flex™ Latex-Free, sterile covers	UA0008*	24 pcs. 11.97 cm narrowing down to 4.57 cm.
Sterile Latex Transducer Covers	UA0005	24 pcs. sterile transducer cover 2.6 x 20 cm.
NeoGuard® SurgiTip Covers, sterile	UA0073	12 pcs., latex-free, 15.2 x 244 cm with SurgiTip. 3 cm, telescopically folded (recommended).

* Must be used together



T-Shaped Intraoperative Transducer 8816

Systems: 2202, 1202-3, 1202-4 and 1202-7
Frequency Range: 10-4 MHz

Carrying Case	KE4300	
Insert for Carrying Case	DS0976	
Transducer Holder	DZ9756	
Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004*	24 pcs. 2.5 m.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009*	24 pcs. 8.89 x 60 cm.
NeoGuard® SurgiTip Covers, sterile	UA0073	12 pcs., latex-free, 15.2 x 244 cm with SurgiTip. 3 cm, telescopically folded (recommended).

* Must be used together

Transducers and Transducer Accessories

Biplane and Endfire Transducer 8818



Systems: 2202, 1202-1 and 1202-3

Frequency Range: 12-4 MHz



UA1329-S

Dummy Channel Bracket	UA1325	Included.
Carrying Case	KE4300	Must be ordered separately.
Civco Multi-Purpose Workstation	UA2020	For transperineal biopsies.
Endfire Biopsy Guides*	UA1327	Reusable. Bore diameter 1.6 mm; angle of insertion 0° to the long axis of the transducer.
	UA1323-S	Single-use, sterile packed. 18 pcs.
Simultaneous Biplane Biopsy Guides*	UA1326	Reusable. Bore diameter 1.6 mm; angle of insertion 19° to the long axis of the transducer.
	UA1322-S	Single-use, sterile packed. 18 pcs.
	UA1322-S14	Sterile packed. 18 pcs. Bore diameter of 2.1 mm; angle of insertion 19° to the long axis of the transducer.
Dual Biopsy Guides*	UA1328	Reusable. Bore diameter 1.6 mm; angle of insertion 0° and 19° to the long axis of the transducer.
	UA1329-S	Single-use, sterile packed. 18 pcs.
Perineal Puncture Guide	UA1324	Perineal puncture attachment. Bore diameter 2.1 mm; angle of insertion 0° to the long axis of the transducer.
Single-Use Cleaning Brush, thin	QZ0040	For needle guide.
Leakage Testing Kit	UA1404	
Magnetic Wheel Mover for 3D Imaging	UA0513	
Mover Adapter Ring	UA0512	
Variable Friction Support Arm	UA0553	For Magnetic Wheel Mover.
Pullback Adapter	UA0516	For Magnetic Wheel Mover.
CIV-Flex™ Latex-Free, sterile covers	UA0008	24 pcs. 11.97 cm narrowing down to 4.57 cm.
Latex Covers	UA0032	24 pcs., non-sterile.
Latex Covers	UA0069	24 pcs., sterile, 3.2 cm tapered to 1.9 x 32 cm.

* All biopsy guides include needle guides

Transducers and Transducer Accessories

Endovaginal Endfire Transducer 8819



Systems: 2202, 1202-1, 1202-2, 1202-3 and 1207-7
Frequency Range: 9-6 MHz



UA1282

Carrying Case	KE4325	
Insert for Carrying Case	DS0537	
Needle Guide	UA1282	Bore diameter 1.8 mm; angle of insertion 0.9° to the long axis of the transducer. Includes 3 single-use cleaning brushes.
Single-Use Cleaning Brush, thin	QZ0040	For needle guide.
CIVCO Sterile Endocavity Needle Guide	610-958	To purchase, please contact your CIVCO supplier.
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
Sterile Latex Transducer Covers	UA0005	24 pcs. sterile covers 2.6 x 30 cm.
CIV-Flex™ Latex-Free, Sterile Covers	UA0008	24 pcs. 11.97 cm narrowing down to 4.57 cm.

Convex Array Transducer 8820e



Systems: 2202, 1202-1, 1202-2, 1202-3 and 1202-7
Frequency Range: 6-2 MHz



UA1250



UA1341 with UA0013

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.
Needle Guide Kit	UA1341	Including: 1 reusable bracket, 5 sterile Ultra-Pro 3™ needle guides with transducer cover and gel.
Sterile Ultra-Pro 3™ Needle Guides	UA0013	24 pcs.
Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009	24 pcs. 8.89 x 60 cm.
CIV-Flex™ Covers	UA0007	24 pcs. sterile covers 13.97 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0070	24 pcs., 14 x 91 cm, telescopically folded
Sterile Latex Covers, intraoperative	UA0072	12 pcs., 20.3 x 244 cm with 9 cm tip, telescopically folded.
CIV-Flex™ Covers, latex-free	UA0076	24 pcs., 14 x 147 cm, telescopically folded.

Transducers and Transducer Accessories

Vascular Transducer 8822



Systems: 2202, 1202-1 and 1202-3
Frequency Range: 9-3.5 MHz

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1239	Bore diameter up to 2.1 mm and angle of insertion 30°, 45° and 60° to the image axis.
Leakage Testing Kit	UA1404	



UA1239

CIV-Flex™ Covers	UA0007	24 pcs. sterile covers 13.97 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0070	24 pcs., 14 x 91 cm, telescopically folded.
NeoGuard® Sterile Covers	UA0074	12 pcs., latex-free, 20.3 x 244 cm with SurgiTip 6 cm, telescopically folded.

Convex Array Transducer 8823



Systems: 2202 and 1202
Frequency Range: 6-1.8 MHz

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.



UA1250

Needle Guide Kit	UA1341	Including: 1 reusable bracket, 5 sterile Ultra-Pro 3™ needle guides with transducer cover and gel.
Sterile Ultra-Pro 3™ Needle Guides	UA0013	24 pcs.



UA1341 with UA0013

Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
CIV-Flex™ sterile, Latex-Free Covers	UA0009	24 pcs. 8.89 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0071	24 pcs., 8.9 x 91.5 cm, telescopically folded.
CIV-Flex™ Cover, latex-free	UA0076	24 pcs., 14 x 147 cm, telescopically folded.

Transducers and Transducer Accessories



I and T Shaped Intraoperative Transducer 8824

Systems: 2202 and 1202-3
Frequency Range: 10-5 MHz

Carrying Case	KP0212	Must be ordered separately.
Insert for Carrying Case	DS0977	
Single-Use Needle Guide Kit UA1337 includes:	UA1337	Box of 9. 2 needle guide insert palettes (3-angle and free-angle respectively), each with 9 different inserts for 10-20 gauge needles.
Leakage Testing Kit	UA1404	
Sterile Cable Covers, latex-free	UA0004*	24 pcs. 2.5 m.
CIV-Flex™ Latex-Free, sterile covers	UA0008*	24 pcs. 11.97 cm narrowing down to 4.57 cm.
NeoGuard® SurgiTip Covers, sterile	UA0073	12 pcs., latex-free, 15.2 x 244 cm with SurgiTip 3 cm, telescopically folded (recommended).



UA1337

* Must be used together



ProART Robotic Transducer 8826

Systems: 2202-7 and 1202-3
Frequency Range: 12-4.3 MHz (depending on system)

Carrying Case	KE4320	Included.
Leakage Testing Kit	UA1404	
Latex-Free Laparoscopic Transducer Cover	UA0067	12 pcs. Telescopically folded (not suitable for robotic surgery).

Transducers and Transducer Accessories



Convex Array Transducer 8830

Systems: 1202
Frequency Range: 6-2 MHz



UA1250



UA1341 with UA0013

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.
Needle Guide Kit	UA1341	Including: 1 reusable bracket, 5 sterile Ultra-Pro 3™ needle guides with transducer cover and gel.
Sterile Ultra-Pro 3™ Needle Guides	UA0013	24 pcs.
Sterile Cable Covers, latex-free	UA0004	24 pcs., 2.5 m.
CIV-Flex™ Sterile Covers	UA0007	24 pcs., 2.5 m.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009	24 pcs., 13.97 x 60 cm.
CIV-Flex™ Sterile Covers, latex-free	UA0070	24 pcs., 8.89 x 60 cm.
Sterile Latex Covers, intraoperative	UA0072	24 pcs., 14 x 91 cm, telescopically folded. 12 pcs., 20.3 x 244 cm with 9 cm tip, telescopically folded.
CIV-Flex™ Covers, latex-free	UA0076	24 pcs., 14 x 147 cm, telescopically folded.

Transducers and Transducer Accessories



3DART Transducer 8838

Systems: 2202 with 3D Professional UA2222 and Single Element Module UA2223, 1202-1 and 1202-3 with 3D Professional UA1222

Frequency Range: 12-4 MHz

Carrying Case	KE4325	Included.
Insert for Carrying Case	DS0980	Included.
Endo Transducer Holder Insert for Pro Focus 2202 UltraView	UA1347	
Transducer Holder for Flex Focus 1202	UA1308	
Robotic Stationary Transducer (RST) TM Arm	UD1309	Designed to hold the transducer during robotic-assisted surgery procedures.
General Transducer Holder	UA1296	For Flex Focus 1202 Mobile Touch Docks UA1810 and UA1814.
Sterile Cable Covers, latex-free	UA0004	24 pcs., 2.5 m.
Sterile Latex Transducer Covers	UA0006	24 pcs. sterile covers 3.5 x 20 cm.
Sterile Latex Transducer Covers	UA0069	24 pcs., sterile, 3.2 cm tapered to 1.9 x 32 cm.

Transducers and Transducer Accessories

Transducer 8848



Systems: 2202, 1202-1 and 1202-3
Frequency Range: 12-4 MHz



UA1232

Carrying Case	KE4326	Must be ordered separately.
Needle Guide	UA1232	Bore diameter variable of 2.2 mm and angle of insertion 0° to the long axis of the transducer.
CIV Flex Brachy Balloons, latex-free	UA0059	10 pcs.
CIVCO Classic Stepper and Micro-Touch® Stand	UA2000	See page 39 for more details.
CIVCO Classic Stepper and Micro-Touch® Stand	UA2001	See page 39 for more details.
EXII Tracked Stepper and Stand	UA2010	See page 39 for more details.
EXII Tracked Stepper and Stand	UA2011	See page 39 for more details.
Magnetic Wheel Mover for 3D Imaging	UA0513	For Magnetic Wheel Mover. For Magnetic Wheel Mover.
Mover Adapter Ring	UA0510	
Variable Friction Support Arm	UA0553	
Pullback Adapter	UA0515	
Sterile Latex Transducer Covers	UA0005	24 pcs. sterile transducer cover 2.6 x 20 cm.
CIV-Flex™ Latex-Free, sterile covers	UA0008	24 pcs. 11.97 cm narrowing down to 4.57 cm.

Craniotomy Transducer 8862



Systems: 2202, 1202-1 and 1202-3
Frequency Range: 10-5 MHz

Carrying Case	KE4300	Must be ordered separately.
Transducer Holder	DZ9756	Included.
Assembly Guide for UA1345	BB2080	Single-Use, Sterile Packed, Needle Guide Kit (box of 9). Includes 1 needle guide and 2 needle guide insert palettes (3-angle and free-angle respectively), each with 9 different inserts for 10-20 gauge needles.
Needle Guide	UA1345	
NeoGuard® SurgiTip Covers, sterile	UA0073	12 pcs., latex-free, 15.2 x 244 cm with SurgiTip 3 cm, telescopically folded (recommended).

Transducers and Transducer Accessories



Burr Hole Transducer 8863

Systems: 1202-3 and 2202-7
Frequency Range: 10-5 MHz

Carrying Case	KE4300	Must be ordered separately.
Transducer Holder	DZ9756	Included.
Assembly guide for UA1346*	BB2080	
Needle Guide	UA1346	Single-Use, Sterile Packed, Needle Guide Kit (box of 9). Includes 1 palette with needle guide and 9 different inserts for 10-20 gauge needles.
Leakage Testing Kit	UA1404	
NeoGuard® SurgiTip Covers, sterile	UA0073	12 pcs., latex-free, 15.2 x 244 cm with SurgiTip 3 cm, telescopically folded (recommended).

*For future release



High Frequency Linear Array Transducer 8870

Systems: 2202-7, 1202-1, 1202-3 and 1202-4
Frequency Range: 18-6 MHz



UA1251

Carrying Case	KE4300	Must be ordered separately.
Needle Guide	UA1251	Bore diameter variable from 0.6-3.4 mm; angle of insertion of 30, 45 or 60° to the long axis of the transducer.
Sterile Cable Covers, latex-free	UA0004	24 pcs. 2.5 m.
CIV-Flex™ Sterile, Latex-Free Covers	UA0009	24 pcs. covers 13.97 x 60 cm.

Transducers and Transducer Accessories

CIVCO AccuCare Stepper Units (Barzell-Whitmore)

For transducer:
8848



CIVCO AccuCare
Micro-Touch® Stepper

CIVCO Classic Stepper and Micro-Touch® stand. Table mount.	UA2000	Including 2 x 5 18G Templates (UA2005), 20 Sterile Drapes (UA2007) and 10 Non-Sterile Brachy Balloons (UA0059).
CIVCO Classic Stepper and Micro-Touch® stand. Table mount.	UA2001	Including 2 x 5 17G Templates (UA2006), 20 Sterile Drapes (UA2007) and 10 Non-Sterile Brachy Balloons (UA0059).
EXII Tracked Stepper and Micro-Touch® stand. Table mount.	UA2010	Including 2 x 5 18G Templates (UA2005), 20 Sterile Drapes (UA2007) and 10 Non-Sterile Brachy Balloons (UA0059).
EXII Tracked Stepper and Micro-Touch® stand. Table mount.	UA2011	Including 2 x 5 17G Templates (UA2006), 20 Sterile Drapes (UA2007) and 10 Non-Sterile Brachy Balloons (UA0059).
CIVCO Multi-Purpose Workstation	UA2020	

Options:

CIV Flex Brachy Balloons, latex-free	UA0059	10 pcs.
Sterile, Single-Use 18G Templates	UA2005	1 box = 5 units. Min. order 2 boxes.
Sterile, Single-Use 17G Templates	UA2006	1 box = 5 units. Min. order 2 boxes.
Sterile Drapes	UA2007	1 box = 20 units.

Index of Reference Numbers

Systems

1202 **1**
1202-1 **1**
1202-2 **1**
1202-3 **1**
1202-4 **1**
1202-7 **1**
2202 **1**
2202-7 **1**
2202-2 **1**
2202-7 **2**
2202-9 **1**
610-958 **31**

A

AN0080 **14**
AN0081 **14**
AN0082 **14**
AN0083 **14**
AN0084 **14**
AN0085 **14**
AN0086 **14**
AN0087 **14**
AN0088 **14**
AN0089 **14**
AN0090 **14**
AO0421 **7**
AO0491 **5, 11**

B

BB1943 **28**
BB2080 **36, 37**

D

DP0925 **7**
DS 0977 **33**
DS0524 **22**
DS0537 **31**
DS0967 **7**
DS0976 **29**
DS0977 **33**
DS0980 **35**
DZ 9756 **36**
DZ9756 **7, 27, 28, 29, 36, 37**
DZ9772 **8**

E

EA4006 **20**
EA4020 **19**
EL4036 **12**
EQ4042 **13**
EQ4072 **13**
EQ4093 **13**
EQ4100 **13**
EQ4102 **8**

EQ4104 **13**
EQ4105 **13**

K

KE1325 **26, 30**
KE4300 **22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 34, 36, 37**
KE4320 **22, 33**
KE4322 **21**
KE4325 **31, 35**
KE4326 **36**
KP0212 **33**

N

NP0959 **8**

Q

QA0185 **22**
QA0228 **5, 7**
QS0011 **20**
QZ 0040 **26**
QZ0040 **22, 25, 30, 31**

S

SC1474 **18**
SC1624 **8**

U

UA0004 **22, 23, 24, 27, 28, 29, 31, 32, 33, 34, 35, 36, 37**
UA0005 **21, 29, 31, 36, 37**
UA0006 **22, 23, 35**
UA0007 **27, 31, 32, 34**
UA0008 **23, 25, 26, 27, 29, 30, 31, 33, 36**
UA0009 **23, 24, 28, 29, 31, 32, 34, 37**
UA0013 **31, 32, 34**
UA0028 **23, 28**
UA0032 **25, 26, 30**
UA0037 **21**
UA0039 **22**
UA0040 **22**
UA0058 **38**
UA0059 **36, 38**
UA0061 **21**
UA0067 **22, 33**
UA0069 **20, 22, 25, 26, 30, 35**
UA0070 **27, 31, 32, 34**
UA0071 **23, 24, 28, 32**
UA0072 **31, 34**
UA0073 **29, 33, 36, 37**
UA0074 **23, 27, 28, 32**
UA0075 **27**
UA0076 **7, 24, 31, 32, 34**
UA0210 **25**
UA0510 **12, 36**
UA0511 **12, 25**
UA0512 **12, 26, 30**

Index of Reference Numbers

UA0513	12, 25, 26, 30, 36	UA1272	25
UA0515	12, 36	UA1282	37
UA0516	12, 26, 30	UA1283	4, 9
UA0553	12, 21, 25, 26, 30, 36	UA1284	10
UA0610	18	UA1284-K	5
UA0650	19, 20	UA1285	5, 10
UA0652	19	UA1294	4, 10
UA0670	12, 21	UA1295	4, 10, 11
UA0671	21	UA1296	5, 35
UA0672	21	UA1298	4
UA0673	21	UA1302	8
UA0674	21	UA1308	35
UA0675	21	UA1322-S	26, 30
UA0676	21	UA1322-S14	26, 30
UA0677	21	UA1323-S	30
UA0678	21	UA1324	26, 30
UA0679	21	UA1326	26, 30
UA0681	19, 21	UA1327	30
UA0682	19, 21	UA1328	30
UA0799	19, 20, 22, 25	UA1329-S	30
UA0878	19	UA1329-S1	30
UA0955	8	UA1335	28
UA0964	5	UA1336	29
UA1166	22	UA1337	33
UA1202	5, 7, 11, 12	UA1341	31, 32, 34
UA1203	4	UA1345	36
UA1204	4, 10	UA1346	37
UA1206	4	UA1347	35
UA1207	4	UA1404	22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33, 37
UA1210	4, 9	UA1414	33
UA1211	4, 9	UA1810	4, 9
UA1214	9	UA1811	4, 9
UA1215	7, 12	UA1814	4, 10
UA1220	5, 10	UA2000	36, 38
UA1222	4	UA2001	36, 38
UA1223	5	UA2005	38
UA1225	5	UA2006	38
UA1226	4	UA2007	38
UA1232	25, 36	UA2010	36, 38
UA1237	4, 10	UA2011	36, 38
UA1239	27, 32	UA2020	26, 30, 38
UA1246	23, 28	UA2204	7
UA1247	5, 9	UA2205	2, 4
UA1249	22	UA2206	2
UA1250	24, 31, 32, 34	UA2207	2
UA1251	37	UA2214	2
UA1256	23	UA2215	2
UA1256-U	23	UA2220	2
UA1257	25	UA2222	2
UA1257-R14	25	UA2222-K	2
UA1257-R17	25	UA2223	2
UA1257R17E	25	UA2223-K	2
UA1257-S14	25	UA2226	2
UA1257S17E	25	UA2226-K	2
UA1258	4	UA2229-K	3
UA1262	27	UA2233-K	3
UA1269	23	UA2235	2
UA1271	7	UA2235-K	2

Index of Reference Numbers

UA2236 **2**
UA2236-K **2**
UA2237-K **2**
UA2238 **2**
UA2238-K **2**
UA2240 **8**
UA2251 **2**
UA2251-K **2**
UA2252 **2**
UA2252-K **2**
UA2410 **3, 7**
UA4100 **3, 8**
UA4102 **5, 13**
UC5296 **25**
UC5296-E **25**
UC5298 **25**
UC5299-E **22**
UD1309 **35**

W

WA0499 **19**
WA0500 **19**
WA0543 **19**
WA0597 **19**
WQ0410 **19**
WQ0739 **18**
WQ0970 **18**
WQ0971 **18**
WQ0974 **18**
WQ0975 **18**
WS3416 **19**

Y

YA0514 **19**
YJ0619 **20**
YJ0621 **19**
YJ0632 **20**
YM1615 **19**

To order, Please contact

MANUFACTURER

BK Medical ApS
Mileparken 34
DK-2730 Herlev
Denmark

Headquarters USA:

B-K Medical Systems, Inc.
8 Centennial Drive
Peabody, MA 001960
USA
Tel.: (1) 978-326-1300
Fax: (1) 978-326-1399

B-K Medical Benelux
Generaal de Wittelaan 19/7
2800 Mechelen
Belgium
Tel.: (32) 1528 0600
Fax: (32) 1529 1025

B-K Medical Medizinische systeme GmbH
Pascalkehe 13
25451 Quickborn
Germany
Tel: (49) 4106 99 55-0
Fax: (49) 4106 99 55-99

B-K Medicale S.r.l.
Via Tulipani, 3/B
20090 Pieve Emanuele MI
Italy
Tel.: (39) 02 90 78 13 47
Fax: (39) 02 90 78 19 05

B-K Medical (UK)
11 Grove Park
Waltham Road, White Waltham
Berkshire SL6 3LW
United Kingdom
Tel.: (44) 1628 825-770
Fax: (44) 1628 826 970

B-K Medical AB
Mänskärsvangen p
141 75 KUNGENS KURVA
Sweden
Tel.: (46) 8-744 02 11
Fax: (46) 8-744 02 12

Your local distributor



Headquarters USA · B-K Medical Systems, Inc. · 8 Centennial Drive, Peabody, MA 01960 · Tel: +1 800 876-7226 / +1 978 326-1300 ■ Headquarters Europe and Rest of World · BK Medical · Mileparken 34, DK-2730 Herlev, Denmark · Tel.: +45 4452 8100
www.bkmed.com

A wholly owned subsidiary of Analogic Corporation