

Transducers and Transducer Accessories



Transducer 8802

Scanners: 2101, 2102 series, 2202 and 2400
Frequency Range: 3- 6 MHz



UA 1250

Carrying Case	KE 4300	Must be ordered separately
Puncture attachment	UA 1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.
Leakage Testing Kit	UA 1404	
Cable covers	UA 0004	24 pcs. sterile cable covers 2.5 m.
CIV-Flex™ covers	UA 0009	24 sterile covers 8.89 x 60 cm.
CIV-Flex™ sterile covers, Latex-free	UA 0071	24 pcs., 8.9 x 91.5 cm, telescopically folded
CIV-Flex™ cover, Latex-free	UA 0076	24 pcs., 14 x 147 cm, telescopically folded



Transducer 8803

Scanners: 1402, 2101, 2102 series, 2202 and 2400
Frequency Range: 3 - 5 MHz



UA 1250

Carrying Case	KE 4300	Must be ordered separately
Puncture attachment	UA 1250	Bore diameter variable from 0.6 to 2.4 mm and angle of insertion 18° to the long axis of the transducer.
Leakage Testing Kit	UA 1404	
Cable covers	UA 0004	24 pcs. sterile cable covers 2.5 m.
CIV-Flex™ covers	UA 0009	24 sterile covers 8.89 x 60 cm.
CIV-Flex™ sterile covers, Latex-free	UA 0071	24 pcs., 8.9 x 91.5 cm, telescopically folded
CIV-Flex™ cover, Latex-free	UA 0076	24 pcs., 14 x 147 cm, telescopically folded